





CEREC...Simplified, Precise & Aesthetic



The CEREC technology provides an ideal way to precisely restore teeth with restorations that fit like gold. It is as simple as capture, design and bond. Patient's love the conservative approach and completing the treatment in one appointment.



With the addition of a digital bite registration, the Biogeneric CEREC software now proposes 95% or more of the final restoration design. This significantly simplifies the design process for even the most computer novice. Restorations can be milled in 5-8 minutes. Glazing and some staining characterization maximizes function and aesthetics.



It is customary to try-in ceramics prior to bonding. Note the accuracy of the dry-fit restoration. This rivals the fit accomplished by the very best labs in the country. Restoration bonded with Excite and Variolink II Transparent. Note the melting of the Triluxe margins. CEREC restorations fit like gold and can be very aesthetic.